ENGLISH

MONTHLY GENERAL CLINICAL CHEMISTRY PROGRAMME: RQ9128

CONFIRMATION OF KIT CHARACTERISTICS AND RECEIPT DATE

Please confirm that the correct number of samples are present and that your samples have the appearance as indicated in the CHARACTERISTICS section below. Please confirm that none of the vials are broken and notify your local Randox representative immediately if there are any discrepancies. Finally, please log on to www.riqas.net to confirm the exact date on which you received this kit.

CHARACTERISTICS

The pack contains 6 vials of lyophilised 5 ml samples. The vials are labelled with the sample number.

PREPARATION

The samples are sealed under vacuum. Open the vial very carefully, avoiding any loss of material and using a **calibrated pipette**, reconstitute with an accurately measured **5 ml volume** of freshly double distilled water at **+20°C to +25°C**. Replace the rubber stopper, close vial and leave to stand for **60 minutes out of bright light** before use. Ensure that the lyophilised contents are completely dissolved by swirling gently. Do not shake the vial. Please analyse the sample within 2 days of reconstitution in the absence of bacterial contamination. The samples should be stored at 2 - 8°C when not in use. The samples should be treated in the same way as patient samples, with the following exceptions.

ANALYSIS OF SAMPLES

Prostatic and Total Acid Phosphatase: The samples should be stabilised by adding 1 drop of acetic acid solution provided to 1 ml of serum.

Bilirubin in the serum is light sensitive so the serum should be **stored in the dark**.

NEFA is stable for 24 hours at +2°C to +8°C following reconstitution.

Bacterial contamination of the reconstituted serum will cause reductions in the stability of many components.

SAFETY

The serum is human based. It has been source tested and found to be negative for HBsAg and antibodies to HIV. For complete protection however, it is recommended that the serum be handled as carefully as patient specimens.

For **IN VITRO** use only.

*** IMPORTANT NOTE:** Results must arrive at **RIQAS** by **17:00 hrs GMT** on the **FINAL DATE**. If the **RECOMMENDED ANALYSIS DATE** gives insufficient time, we suggest that the sample is analysed earlier to ensure you meet the deadline. If you are faxing results, please transmit **3 working days before the FINAL DATE**. Late results will not be accepted after the final date for the next sample.

MONTHLY GENERAL CLINICAL CHEMISTRY PROGRAMME / PROGRAMME CHIMIE CLINIQUE GENERAL MENSUEL /
 PROGRAMMA MENSILE DI CHIMICA CLINICA GENERALE / PROGRAMA MENSUAL DE QUÍMICA CLÍNICA GENERAL /
 PROGRAMA DE QUÍMICA CLÍNICA GERAL MENSAL / MÁNATLIGT ALLMÁNT KLINISKT KEMISKT PROGRAM /
 按月临床化学分析程序 / CHƯƠNG TRÌNH NGOẠI KIỂM SINH HÓA

RETURN OF RESULTS / DATES de RETOUR / INVIO DE RISULTATI / ENVIO DE RESULTADOS /
 PLAN FOR INSENDELSE AF / RESULTATRETUR / 回复结果 / GỬI TRẢ KẾT QUẢ /

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