ENGLISH

GLYCATED HAEMOGLOBIN (HbA1c) PROGRAMME: RQ9129

CONFIRMATION OF KIT CHARACTERISTICS AND RECEIPT DATE

Please confirm that the correct number of samples are present and that your samples have the appearance as indicated in the CHARACTERISTICS section below. Please confirm that none of the vials are broken and notify your local Randox representative immediately if there are any discrepancies. Finally, please log on to www.riqas.net to confirm the exact date on which you received this kit.

CHARACTERISTICS

The pack contains 6 vials of lyophilised 0.5 ml samples. The vials are labelled with the sample number.

PREPARATION

The samples are sealed under vacuum. Open the vial very carefully, avoiding any loss of material and reconstitute with an accurately measured 0.5ml volume of freshly distilled water at +20°C to +25°C. Replace the sample bottle cap, swirl the bottle several times and leave to stand at room temperature for 15 minutes. After 15 minutes, coat all surfaces of the bottle by rotating and inverting the bottle. Continue mixing until the solution is homogeneous and all lyophilised material is reconstituted. Do not shake the vial. The reconstituted sample should be analysed within 2 days. Do not freeze the reconstituted samples. The samples should be stored at 2 to 8°C when not in use.

NB. The reconstituted samples should be treated the same as patient samples. Depending on the method of analysis, they may require pre-treatment (either on the analyser, or manually) in order to assay for HbA1c.

SAFETY

Potentially Biohazardous Material. The samples contain human blood. Human source material from which this product has been derived has been tested at donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B surface Antigen (HbsAg), and Hepatitis C Virus (HCV) antibody and found to be NON-REACTIVE. FDA approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, this material and all patient samples should be handled as though capable of transmitting infectious diseases and disposed of accordingly.

For **IN VITRO** use only, do not pipette by mouth, exercise the normal precautions required for handling laboratory reagents.

*** IMPORTANT NOTE:** Results must arrive at **RIQAS** by **17:00 hrs GMT** on the **FINAL DATE**. If the **RECOMMENDED ANALYSIS DATE** gives insufficient time, we suggest that the sample is analysed earlier to ensure you meet the deadline. If you are faxing results, please transmit **3 working days before the FINAL DATE**. Late results will not be accepted after the final date for the next sample.

GLYCATED HAEMOGLOBIN (HbA1c) PROGRAMME / PROGRAMME HEMOGLOBINE GLYQUEE (HbA1c) /
 PROGRAMMA EMOGLOBINA GLICATA (HbA1c) / PROGRAMA DE HEMOGLOBINA GLICOSILADA (HbA1c) /
 HEMOGLOBINA GLICADA (HbA1c) PROGRAMA / GLYCATED HAEMOGLOBIN (HbA1c) PROGRAM /
 PROGRAM DO OZNACZANIA HEMOGLOBINY GLIKOWANEJ (HbA1c) / CHƯƠNG TRÌNH NGOẠI KIỂM (HbA1c)

RETURN OF RESULTS / DATES de RETOUR / INVIO DEI RISULTATI / ENVIO DE RESULTADOS
 PLAN FOR INSENDELSE AF / ZWROTU WYNIKÓW / GỬI TRẢ KẾT QUẢ

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