ENGLISH
BLOOD GAS PROGRAMME: RQ9134

CONFIRMATION OF KIT CHARACTERISTICS AND RECEIPT DATE

Please confirm that the correct number of samples are present and that your samples have the appearance as indicated in the CHARACTERISTICS section below. Please confirm that none of the vials are broken and notify your local Randox representative immediately if there are any discrepancies. Finally, please log on to www.riqas.net to confirm the exact date on which you received this kit.

CHARACTERISTICS

The pack contains 6 vials of liquid material (6 x 1.8 ml). The vials are labelled with the sample number.

STORAGE

The samples should be stored at 2-8°C. Avoid exposure to freezing and temperatures greater than 30°C.

PREPARATION

The samples should be brought to a temperature of 20 – 23°C before use. Allow at least 4 hours for ampoules to equilibrate to this temperature prior to testing. (Note: pO₂ values will vary inversely by about 1% per degree Celsius that the temperature of the ampoules varies). For pH/blood gas values, the samples should be analyzed within 1 minute of opening. For electrolyte measurements the samples should be analyzed within 1 hour after opening. Follow the procedures listed below:

1. Before use, hold the ampoule at the top and bottom (with forefinger and thumb) and shake 15-20 times (about 10 seconds) to mix the solution. Tap the ampoule to restore the liquid to the bottom of the ampoule.
2. Open the ampoule by snapping off the tip at the score. Use gauze, tissue, gloves, or an appropriate ampoule opener to protect fingers from cuts.
3. Immediately introduce the liquid from the ampoule to the analyser. Follow the manufacturer's instructions for sampling a control material. Depending on the sampling procedure chosen, the following instructions apply:
 - a. Direct Aspiration: Sample directly from the ampoule.
 - b. Syringe Transfer:
 - i. Use a clean, gas-tight syringe attached to a clean, blunt syringe needle (if available)
 - ii. Prime the syringe by slowly aspirating a small amount (0.2-0.3 ml) of solution from the ampoule.
 - iii. Discard this liquid, leaving the dead space of the syringe filled with the sample.
 - iv. Aspirate the sample from the ampoule into the primed syringe. Be careful that air is not drawn in with the liquid. Expel 1 to 2 drops, detach the needle and immediately inject the sample into the analyzer sample port.
 - c. Ampoule Injector/Dispenser: Assemble and fill the ampoule injector following the manufacturer's instructions. Expel 1 or 2 drops to rinse the outlet tip and inject the sample into the analyzer sample port.
 - d. Capillary Mode:
 - i. Install the appropriate adapter for micro sampling onto the instrument.
 - ii. Sample the contents of the ampoule following the recommendations of the instrument manufacturer. Be certain to keep the sampling tip of the adapter below the surface of the liquid during aspiration.

The samples should be treated in the same way as patient samples. If samples are being run on any of the following Roche instruments **cobas b 121**, **cobas b 121 BGE**, **cobas b 221**, **cobas b 123** the samples should be measured as **normal blood samples**. **Riqas Samples should be run on the GEM Premier 4000 in "Proficiency Mode" only**. Those using a Siemens Rapid Blood Gas analyser should measure samples as 'patient samples' and run in syringe mode for arterial blood.

SAFETY

These samples contain no human or biological materials. For IN VITRO use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.

*** IMPORTANT NOTE:** Results must arrive at **RIQAS** by **17:00 hrs GMT** on the **FINAL DATE**. If the **RECOMMENDED ANALYSIS DATE** gives insufficient time, we suggest that the sample is analysed earlier to ensure you meet the deadline. If you are faxing results, please transmit **3 working days before the FINAL DATE**. Late results will not be accepted after the final date for the next sample.

BLOOD GASES PROGRAMME / PROGRAMME GAZ DU SANG / PROGRAMMA GAS SANGUIGNI / PROGRAMA GAS ARTERIALES / KRAUJO DUJŲ PROGRAMA / โปรแกรม BLOOD GAS / PROGRAM GAZOMETRIA / 血气质控程序 / CHƯƠNG TRÌNH NGOẠI KIỂM KHÍ MÁU

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