ENGLISH
LIQUID CARDIAC PROGRAMME: RQ9136

CONFIRMATION OF KIT CHARACTERISTICS AND RECEIPT DATE

Please confirm that the correct number of samples are present and that your samples have the appearance as indicated in the CHARACTERISTICS section below. Please confirm that none of the vials are broken and notify your local Randox representative immediately if there are any discrepancies. Finally, please log on to www.riqas.net to confirm the exact date on which you received this kit.

CHARACTERISTICS

The pack contains 6 vials of liquid sample (6 x 3 ml). The vials are labelled with the sample number.

PREPARATION

Ensure that samples are homogeneous by swirling gently. Do not shake the vials. Once opened the samples should be analysed within 2 days. The samples should be treated in the same way as patient samples. The samples should be stored at 2 - 8°C when not in use.

SAFETY

Warning: Potentially Biohazardous Material

Human source material from which this product has been derived has been tested at donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B surface Antigen (HbsAg), and Hepatitis C Virus (HCV) antibody and found to be NON-REACTIVE. FDA approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, this material and all patient samples should be handled as though capable of transmitting infectious diseases and disposed of accordingly.

The samples contain < 0.1% w/v Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention. Sodium Azide reacts with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide.

For IN VITRO use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.

*** IMPORTANT NOTE:** Results must arrive at **RIQAS** by **17:00 hrs GMT** on the **FINAL DATE**. If the **RECOMMENDED ANALYSIS DATE** gives insufficient time, we suggest that the sample is analysed earlier to ensure you meet the deadline. If you are faxing results, please transmit **3 working days before the FINAL DATE**. Late results will not be accepted after the final date for the next sample.

LIQUID CARDIAC PROGRAMME / PROGRAMME MARQUEURS CARDIAQUES LIQUIDE / PROGRAMMA CARDIACO LIQUIDO /
 PROGRAMA CARDIACO LÍQUIDO / PROGRAMA CARDIACO LÍQUIDO / ŠIRDIES ŽYMENŲ SKYSTOS KONTROLINĖS MEDŽIAGOS PROGRAMA /
 โปรแกรม LIQUID CARDIAC / PROGRAM PŁYNNNA KARDIOLOGIA / 液体心肌质控程序 / CHƯƠNG TRÌNH NGOẠI KIỂM TIM MẠCH

RETURN OF RESULTS / RETOUR DES RESULTATS / INVIO DEI RISULTATI / ENVIO DE RESULTADOS /
 REZULTATŲ GRAŽINIMAS / ตารางการส่งผล RIQAS กลับ / PRZESYŁANIE WYNIKÓW / 回复结果 / GỬI TRẢ KẾT QUẢ

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