## ENGLISH

### IMMUNOASSAY SPECIALITY 2: RQ9142

#### CONFIRMATION OF KIT CHARACTERISTICS AND RECEIPT DATE

Please confirm that the correct number of samples are present and that your samples have the appearance as indicated in the CHARACTERISTICS section below. Please confirm that none of the vials are broken and notify your local Randox representative immediately if there are any discrepancies. Finally, please log on to [www.riqas.net](http://www.riqas.net) to confirm the exact date on which you received this kit.

#### CHARACTERISTICS

The pack contains 6 vials of lyophilised material (6 x 1 ml).

#### PREPARATION

The samples are sealed under vacuum. Open the vial very carefully, avoiding any loss of material and using a **calibrated pipette** reconstitute with an accurately measured **1 ml volume** of freshly double distilled water at **+20°C to +25°C**. Replace the rubber stopper and ensure that samples are dissolved completely by swirling gently (ideally place on a roller for half an hour prior to analysis). Do not shake the vials. Once reconstituted Renin is stable for 5 days when stored at +2 to +8 °C, Prolactin is stable for 1 day and Gastrin and Calcitonin are stable for 8 hours at +2 - 8°C.

The samples should be treated in the same way as patient samples. The samples should be stored at 2 to 8°C when not in use.

#### SAFETY

##### Warning: Potentially Biohazardous Material

Human source material from which this product has been derived has been tested at donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B surface Antigen (HBsAg), and Hepatitis C Virus (HCV) antibody and found to be NON-REACTIVE. FDA approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, this material and all patient samples should be handled as though capable of transmitting infectious diseases and disposed of accordingly.

For IN VITRO use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.

**\* IMPORTANT NOTE:** Results must arrive at **RIQAS** by **17:00 HRS GMT** on the **FINAL DATE**. If the **RECOMMENDED ANALYSIS DATE** gives insufficient time, we suggest that the sample is analysed earlier to ensure you meet the deadline. If you are faxing results, please transmit **3 working days before the FINAL DATE**. **Late results will not be accepted after the final date for the next sample.**

IMMUNOASSAY SPECIALITY 2 PROGRAMME / PROGRAMA DE INMUNOENSAYO ESPECIAL 2/ PROGRAMA DE IMMUNOENSAYO ESPECIALIDAD 2 /  
 PROGRAMMA IMMUNOASSAY SPECIALITY 2 / IMMUNOLOGIA - OZNACZENIA SPECJALISTYCZNE 2 / SPECIALI IMUNOLOGINĚ PROGRAMA 2 /  
 CHƯƠNG TRÌNH NGOẠI KIỂM MIỄN DỊCH ĐẶC BIỆT 2

RETURN OF RESULTS / RENDU DES RESULTATS / DEVOLUCIÓN DE LOS RESULTADOS / ENVIAR RESULTADOS  
 PRZESYŁANIE WYNIKÓW / REZULTATŲ GRAŽINIMAS / GŪI TRĀ KĒT QUĀ

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