**ENGLISH**

**SEROLOGY (ToRCH): RQ9152D**

**CONFIRMATION OF KIT CHARACTERISTICS AND RECEIPT DATE**

Please confirm that the correct number of samples are present and that your samples have the appearance as indicated in the CHARACTERISTICS section below. Please confirm that none of the vials are broken and notify your local Randox representative immediately if there are any discrepancies. Finally, please log on to [www.rigas.net](http://www.rigas.net) to confirm the exact date on which you received this kit.

**CHARACTERISTICS**

20 samples are supplied per cycle and are shipped in packs of 5. The pack contains 5 vials of liquid material (5 x 1 ml). This pack contains samples 16 - 20. The vials are labelled with the sample number.

**PREPARATION AND ANALYSIS OF SAMPLES**

Ensure the samples are homogeneous by swirling gently, do not shake the vials. Centrifuge the samples before analysis at 5000 RPM for 1 minute. If there is insufficient liquid sample remaining, please centrifuge the sample at 15,000 RPM for 5 minutes. To avoid evaporation, close vial immediately after use and store at +2 to +8°C. All samples should be run together and results returned to RIQAS by the date below. The samples should be treated in the same way as patient samples. Once opened please analyse the samples within 2 days.

**SAFETY**

**Warning: Potentially Biohazardous Material. Biological Source Material. Treat as potentially infectious.**

Human source material from which this product has been derived has been tested at donor level for the presence of Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B surface Antigen (HbsAg), and Hepatitis C Virus (HCV) antibody using FDA approved methods and has found to be non-reactive. **However, since no method can offer complete assurance as to the absence of infectious agents, this material should be handled as though capable of transmitting infectious diseases and disposed of accordingly.**

For **IN VITRO** use only.

**\* IMPORTANT NOTE:** Results must arrive at RIQAS by **17:00 HRS GMT** on the **FINAL DATE**. If the **RECOMMENDED ANALYSIS DATE** gives insufficient time, we suggest that the sample is analysed earlier to ensure you meet the deadline. If you are faxing results, please transmit 3 working days before the **FINAL DATE**. Late results will not be accepted after the final date for the next sample.

SEROLOGY (ToRCH) / SEROLOGIE (ToRCH) / SIEROLOGIA (ToRCH) /  
 SEROLOGIA (ToRCH) / 血清 (TORCH) / CHƯƠNG TRÌNH NGOẠI KIỂM ToRCH

RETURN OF RESULTS / RETOUR DES RESULTATS / INVIO DEI RISULTATI / ENVIO DE RESULTADOS /  
 PRZESYŁANIE WYNIKÓW / 回复结果 / GỬI TRẢ KẾT QUẢ

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