English

**ESR PROGRAMME: RQ9163**

**CONFIRMATION OF KIT CHARACTERISTICS AND RECEIPT DATE**

Please confirm that the correct number of samples are present in your kit and that your samples have the appearance as indicated in the CHARACTERISTICS section below. Please confirm that none of the primary tubes are damaged and notify your local Randox representative immediately if there are any discrepancies. Finally, please log on to [www.riqas.net](http://www.riqas.net) to confirm the exact date on which you received this kit.

**IMPORTANT NOTICE**

This ESR control material is not suitable for all ESR analysers. It is not compatible with instruments that use Capillary Photometry Technology that require latex calibrators such as the Alifax. The ESR control material is not suitable for use on the Alcor iSED instrument.

If the tube provided is unsuitable for your instrument, please transfer the sample into your normal testing tube.

**CHARACTERISTICS**

Each pack contains 2 x 4.5 ml of liquid material in primary tubes. The control material is an in vitro reagent composed of mammalian erythrocytes suspended in a plasma like fluid with preservatives. Each primary tube is labeled with the sample number.

**PREPARATION**

**Caution: It is critically important to mix the control material thoroughly at all mixing steps. DO NOT mix on a mechanical mixer.**

**Preparing the control sample:**

- 1) Remove the sample tube to be analysed from the refrigerator and allow the sample to warm to **20 - 25°C** for 15 minutes **before** mixing.

**Mixing the control sample:**

- 2) Mix the sample by holding the tube horizontally between the palms of the hands. Do not pre-mix on a mechanical mixer. Incomplete mixing of the tube prior to analysis invalidates the sample withdrawn and any remaining material in the tube.
  - a) Roll the tube back and forth for 30-60 seconds, occasionally inverting the tube. Mix vigorously but **DO NOT SHAKE**.
  - b) Continue to mix in this manner until the red cells are completely suspended.
  - c) Gently invert the tube 10 times immediately before analysing the sample.
  - d) Analyse the sample immediately after mixing.

**Analysing the control sample:**

- 3) For **automated** methods, **do not** remove the diluent from test reservoirs before using this control.
- 4) For **manual** methods, if you normally dilute patient samples, also dilute this control.

**After sampling:**

- 5) If the tube has been opened for sampling, carefully wipe the rim of the tube and cap with a lint-free tissue and replace the cap tightly.
- 6) Return the tube to the refrigerator within 30 minutes of use.

The samples should be processed and analysed in the same way as patient samples.

**STORAGE AND STABILITY**

Samples should always be stored upright at 2 - 8°C when not in use. **Do not allow tubes to overheat or freeze.** Once opened please analyse the samples within 2 days

**SAFETY**

Potentially hazardous biomaterial. Human Blood components were not used in the manufacture of this control. However, this control does contain components from non-human sources and may transmit infectious disease. When handling or disposing of this product, follow the same precautions that you would for patient samples. Wear protective laboratory gear when handling this control.

For **IN VITRO** use only.

**\* IMPORTANT NOTE:** Results must arrive at RIQAS by **17:00 HRS GMT** on the **FINAL DATE**. If the **RECOMMENDED ANALYSIS DATE** gives insufficient time, we suggest that the sample is analysed earlier to ensure you meet the deadline. If you are faxing results, please transmit 3 working days before the **FINAL DATE**. Late results will not be accepted after the final date for the next sample.

ESR PROGRAMME / PROGRAMME VS/PROGRAMMA VS/ PROGRAMA ESR PROGRAM DE VS/ENG PROGRAM/ PROGRAM OB/ESR PROGRAM		
RETURN OF RESULTS / DATES de RETOUR / INVIO DEI RISULTATI / ENVIO DE RESULTADOS PLAN FOR INDSSENDSE AF/ZWROTU WYNIKÓW		
CYCLE 5A / CICLO 5A / RESULTATER 5A / RESULTATRETUR 5A / CYKL 5A		
SAMPLE NO / ECHANTILLON NO. / N. CAMPIONE / N° DE MUESTRA / AMOSTRA N°. / PRØVENR. / AP. ΔΕΙΓΜΑΤΟΣ / PROV NR. / NR PRÓBK	RECOMMENDED ANALYSIS DATE/ DATE RECOMMANDE POUR L' ANALYSE/ DATA DI ANALISI RACCOMANDATA/ FECHA DE ANALISIS RECOMENDADA/ DATA RECOMENDADA PARA ANÁLISE/ ANBEFALET ANALYSEDATO/ REKOMMENDERAT ANALYSDATUM / ZALECANA DATA WYKONANIA OZNACZENIA	* FINAL DATE/ *DATE LIMITE/ * DATA FINALE/ *FECHA FINAL/ * DATA FINAL/ * SLUTDATO/ * SLUTDATUM/ * DATA KONCOWA
1	18.03.19	25.03.19
2	18.03.19	25.03.19

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