English

BLOOD TYPING PROGRAMME: RQ9206

CONFIRMATION OF KIT CHARACTERISTICS AND RECEIPT DATE

Please confirm that the correct number of samples are present and that your samples have the appearance as indicated in the CHARACTERISTICS section below. Please confirm that none of the vials are broken and notify your local Randox representative immediately if there are any discrepancies. Finally, please log on to www.riqas.net to confirm the exact date on which you received this kit.

CHARACTERISTICS

The kit contains 3 vials. The vials are labelled with the sample number. There will be 4 kits in total throughout the cycle, each sent quarterly.

PREPARATION/STORAGE/STABILITY OF SAMPLES

Samples are ready to use and should be treated as you would treat normal patient/donor samples, while following your laboratory's standard operating procedures and the IFUs of your normal reagents and instrumentation.

Sample 5 = 15% ($\pm 2\%$) Simulated Whole Blood 4 mL – This should be tested for **forward and reverse ABO typing, Rh (D) Typing, antibody screening and identification, and crossmatch with the donor sample.**

Sample 6 = 15% ($\pm 2\%$) Simulated Whole Blood 4 mL – This should be tested for **forward and reverse ABO typing, Rh (D) Typing, antibody screening and identification, and crossmatch with the donor sample.**

Donor = 25% ($\pm 2\%$) RBC Suspension 2 mL – This is the donor sample and should be tested for **compatibility with samples 1 AND 2.** Once opened, the samples should be analysed within 2 days.

SAFETY

Contains human red cells that have been tested for HBsAg, STS, antibodies to HCV, HBc, HIV 1 and 2, West Nile virus, and NAT test for HIV, HBV, and HCV and found acceptable. Each donor is tested for antibodies to T. cruzi once and found acceptable. Contains bovine material obtained from a USDA approved source free of Transmissible Spongiform Encephalopathies (TSEs). However, it is recommended that this product be handled with the same precautions used for patient specimens.

For **IN VITRO** use only.

*** IMPORTANT NOTE:** Results must arrive at RIQAS by **17:00 HRS GMT** on the **FINAL DATE**. If the **RECOMMENDED ANALYSIS DATE** gives insufficient time, we suggest that the sample is analysed earlier to ensure you meet the deadline. Late results will not be accepted after the final date for the next sample.

BLOOD TYPING PROGRAMME // PROGRAMMA GRUPPO SANGUIGNO / PROGRAMA DE TIPIFICACIÓN SANGUÍNEA / / 血型计划 / CHƯƠNG TRÌNH XÁC ĐỊNH NHÓM MÁU / / KAN TÍPI PROGRAMI / BLOOD TYPING PROGRAMME

RETURN OF RESULTS / DATES DE RETOUR / INVIO DEI RISULTATI / ENVIO DE RESULTADOS
 PLAN FOR INSENDELSE AF / CYKEL FÖR / ΚΥΚΛΟΣ ΕΠΙΣΤΡΟΦΗΣ / SONUÇ GÖNDERİM TARİHLERİ / ODOSLANIE VÝSLEDKOV

CYCLE 2C / CICLO 2C / RESULTATER 2C / RESULTATRETUR 2C /
 CYKL ZWROTU WYNIKÓW 2C / 返回结果周期 2C / DÖNEM 2C / CYKLUS 2C

SAMPLE NO / ECHANTILLON NO. /
 N. CAMPIONE / N° DE MUESTRA /
 AMOSTRA N°. / PRØVENR. /
 ΑΡ. ΔΕΙΓΜΑΤΟΣ / PROV NR. /
 NR PRÓBKÍ / 样品号 / ตัวอย่างตรวจที่
 / ÖRNEK NO / Č.VZORKY

RECOMMENDED ANALYSIS DATE/ DATE RECOMMANDE POUR L' ANALYSE/
 DATA DI ANALISI RACCOMANDATA/ FECHA DE ANALISIS RECOMENDADA/
 DATA RECOMENDADA PARA ANÁLISE/ ANBEFALET ANALYSEDATO/
 REKOMMENDERAT ANALYSDATUM/ ΣΥΝΙΣΤΩΜΕΝΗ ΗΜΕΡΟΜΗΝΙΑ ΑΝΑΛΥΣΗΣ
 ZALECANA DATA WYKONANIA OZNACZENIA/ 推荐分析日期 /
 วันที่แนะนำให้ทำการวิเคราะห์ / ÖNERİLEN ANALİZ TARİHİ/
 ODPORÚČANÝ DÁTUM ANALÝZY

* FINAL DATE/ *DATE LIMITE/
 * DATA FINALE/ * FECHA LÍMITE /
 * DATA FINAL/ * SLUTDATO/
 * SLUTDATUM/ *ΤΕΛΙΚΗ ΗΜΕΡΟΜΗΝΙΑ/
 * DATA KOŃCOWA/ *最终日期/
 * วันสุดท้ายของการส่งผลกลับ /
 * SON SONUÇ GÖNDERİM TARİHİ/
 * FINÁLNY DÁTUM

5	20.07.26	27.07.26
6	20.07.26	27.07.26

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